

eas[®]ret

577nm Fiber Technology Laser



Peripheral and Macular Photocoagulation



easyret®


Easyret® is a fully integrated 577nm yellow photocoagulator based on a technological breakthrough: fiber laser technology. Available with Haag Streit or Zeiss type slit lamps, it offers a large choice of treatment settings well adapted to the treatment of macular and peripheral retinal pathologies.

Easyret®: Yellow, MultiSpot and SubLiminal™ Modes

● Yellow Laser - 577nm Wavelength:

Presented as the most versatile wavelength in the scientific literature, the 577nm wavelength offers the following benefits:

- Excellent combined absorption by both melanin and oxyhemoglobin (peak absorption of oxyhemoglobin) [1,2]
- Very little absorption by macular xanthophyll pigments [1,2]
- Excellent penetration through cataracts and hazy media [1,2]



● MultiSpot Mode:

Characterized by the use of short pulse durations from 10 to 20 ms, the MultiSpot treatment mode offers many advantages over classical treatments:

- Less heat diffusion to the retina and choroid, less damage to the retinal nerve fiber layer [3,4]
- Comfortable treatment better tolerated by patients [5]
- Treatment time reduction (full PRP in 1 session) [6]

The MultiSpot treatment mode can be delivered through 5 customizable patterns for better adaptation to the treatment site.

Single spot - Squares - Circles - Triple arcs - Macular grid

PATTERN SELECTION TYPE	DELIVERED LASER SPOTS
SINGLE SPOT	
SQUARES	
CIRCLES	
TRIPLE ARCS	
MACULAR GRID	

● SubLiminal™ Mode:

Composed of a train of extremely short microsecond pulses, this subthreshold treatment mode (non-visible laser impacts) is a tissue sparing treatment mode avoiding scarring [7,8] while treating Diabetic Macular Edema [7] and Central Serous Chorioretinopathy [8].

The SubLiminal™ treatment mode can be delivered through 3 customizable patterns for better adaptation to the treatment site.

PATTERN SELECTION TYPE	DELIVERED LASER SPOTS
SINGLE SPOT	
SQUARES	
CUSTOMIZABLE MACULAR GRID	

Easyret®: Fully Integrated Design

Easyret® offers a fully integrated design in which the laser and the slit lamp are optimally integrated for better ergonomics and ease of use. It is available with two types of slit lamps to adapt to the operator's working habits.

Both versions feature:

- 1 - an integrated laser adapter featuring a continuously variable parfocal zoom.
- 2 - a large touch screen interface to monitor the treatment settings.
- 3 - a click wheel to control the patterns settings
- 4 - an intelligent footswitch to control the laser settings.

Haag Streit Type

Zeiss Type



Easyret®: Enhanced Software User Interface



● 3 Treatment Modes / 3 Dedicated Targets:

Easyret® provides an intuitive and versatile software user interface simplifying the use of the Single Spot, MultiSpot and SubLiminal™ treatment modes. Built in a clinically oriented manner, Easyret® offers 3 different types of visible targets (aiming beam) facilitating the implementation of the laser spots with each treatment mode.



Treatment Report: ●

After treatment, a detailed report can be generated in PDF format. It can be printed and / or saved on a dedicated USB key.





A WORLD FIRST TO MARKET IN PHOTOCOAGULATION: FIBER LASER CAVITY

Easyret®: Technology


● Fiber Laser Technology:

Stemming from the ELBA™ technology, developed and successfully marketed by Quantel Laser for various applications, this new generation of laser cavity provides unique advantages:

- An excellent beam quality ensuring a homogeneous laser spot profile (top hat)
- The emission of pure 577nm yellow wavelength
- An extended lifetime thanks to a simple, compact and reliable technology

The fiber laser technology is a variation of the standard solid-state laser technology.


In fiber lasers, the lasing medium is composed of an optical fiber doped with rare earth elements and optically pumped by diodes.



● Resume® Technology:

Easyret® features the proprietary Resume® function offering more flexibility to the operator in the implementation of the MultiSpot and the SubLiminal™ treatment modes.


- In MultiSpot mode, the pattern delivery can be paused and resumed (the previous shots are remembered)
- In SubLiminal™ mode, the treatment is combined with the pattern scan mode and delivered semi-automatically in several steps



● SubLiminal™ Technology:

In addition to SingleSpot and MultiSpot delivery modes, Easyret® features the SubLiminal™ technology.

The use of this subthreshold treatment mode converts each laser shot into a “pulse envelope” composed of a customizable train of short pulses, allowing the operator to fully adjust the pulse duration (On Time) and interval (Off Time). This finely-tuned control of the laser treatment settings ensures a precise management of the thermal effect on the targeted tissues.



eas^yret[®]

TECHNICAL SPECIFICATIONS



EASYRET SPECIFICATIONS

Laser source:	fiber laser technology
Wavelength:	yellow 577nm
Power at tissue up to:	2000 mW
Pulse duration:	10 ms to continuous
Single spot modes:	single, repeat, painting, continuous
SubLiminal™ mode:	train of microsecond pulses adjustable duty cycle: 5% to 100%
Resume [®] function:	available in Multispot and SubLiminal™ modes
Pattern:	
MultiSpot mode:	single spot, squares, circles, triple arc, macular grid
SubLiminal™ mode:	single spot, squares, customizable macular grid
Spot size:	
Single spot:	continuously variable from 50 µm to 400 µm
Pattern:	continuously variable from 100 µm to 400 µm
Integrated slit lamps:	
Haag Streit type:	Quantel Medical (CSO 9900 5x)
Zeiss type:	Quantel Medical (CSO 9800 5x)
Aiming beam:	635 - 650nm
Size:	174.2 (H) x 97 (W) x 72 (D) cm 68.58" (H) x 38.19" (W) x 28.35" (D)
Weight:	60 kg - 132 lbs
Cooling:	by Peltier effect
Power requirements:	100 to 240 VAC, 350 VA, 50/60 Hz

OPTIONAL FEATURES


Second laser port	
Laser indirect ophthalmoscopes:	Heine Omega 500 or Keeler Vantage Plus

Specifications are subject to change without notice.
©2017, Quantel Medical, Easyret and Resume Function are registered trademarks of Quantel Medical.
Elba is a trademark of Quantel. All rights reserved.

- 1- Vogel M, Schäfer FP, Stuke M, Müller K, Theuring S, Morawietz A. Animal, experiments for the determination of an optimal wavelength for retinal coagulations. Graefes Arch Clin Exp Ophthalmol. 1989;227:277-280.
- 2- Mainster MA. Wavelength selection in macular photocoagulation. Tissue optics, thermal effects, and laser systems. Ophthalmology. 1986;93:952-958.
- 3- Jain A, Blumenkranz MS, Paulus Y et al. Effect of pulse duration on size and character of the lesion in retinal photocoagulation. Arch Ophthalmol. 2008; 126:78-85.
- 4- Yi-Ryeung Park, Donghyun Jee. Changes in Peripapillary Retinal Nerve Fiber Layer Thickness after Pattern Scanning Laser Photocoagulation in Patients with Diabetic Retinopathy. Korean J Ophthalmol 2014;28(3):220-225.
- 5- Hussainy S AI, Dodson PM and Gibson JM Pain response and follow-up of patients undergoing panretinal laser photocoagulation with reduced exposure times. Eye (2008) 22, 96-99
- 6- Muqit MM, Marcellino GR, Henson DB et al. Single-Session vs Multiple-Session Pattern Scanning Laser Panretinal Photocoagulation in Proliferative Diabetic. Arch ophthalmol, 2010, 128: 525-533
- 7- Yoon Hyung Kwon, Dong Kyu Lee, Oh Woong Kwon The short-term efficacy of subthreshold micropulse yellow(577-nm) laser photocoagulation for diabetic macular edema. Korean J Ophthalmol 2014;28(5):379-385
- 8- Scholz P, Ersoy L, Boon CJF, Fauser S Subthreshold Micropulse Laser(577 nm). Treatment in Chronic Central Serous Choriorretinopathy. Ophthalmologica 2015 DOI: 10.1159/000439600

www.quantel-medical.com

A product by **Quantel Medical, France**



CE
0459
ISO 9001- ISO 13485

Headquarters
Quantel Medical
11, rue du Bois Joli - CS40015
63808 Cournon d'Auvergne Cedex - FRANCE
Tel: +33 (0)4 73 745 745
Fax: +33 (0)4 73 745 700
E-mail: contact@quantel-medical.fr

North America
Quantel USA
49 Willow Peak Dr.
Bozeman, MT 59718 – USA
Tel: +1 877 782 6835
Fax: +1 406 522 2005
E-mail: info@quantelmedical.com

Representative Offices
Thailand, Chiang Mai
Brazil, Rio De Janeiro

Quantel
medical