

## 51-935 Artificial Anterior Chamber

The CORONET Artificial Anterior Chamber has been designed to be used with any of the CORONET standard hand held trephines and is unique in concept in that one chamber fits all sizes of donor graft - ONE SIZE FITS ALL.

Suitable for use with both manual dissection and femto-second lasers.

- Supplied one per box
- Supplied sterile for single use only
- Sterilisation by Ethylene Oxide (EO)
- Declared 3-year shelf life

- The chamber is suitable for any size donor graft tissue from 14.5mm through to 18mm
- Two tissue retaining heads are supplied with the chamber:

**Small** to accommodate graft tissue  
14.50mm - 16.00mm

**Large** to accommodate graft tissue  
16.00mm - 18.00mm

