



Product Catalog



Single-Use Ophthalmic Surgical Products



How to Order

Toll Free Phone: 800 782-6534
Phone: 909 481-0011
Fax: 909 481-4481

Please contact us via phone, fax, email or mail and you will receive prompt assistance from a qualified customer service representative.

Order online: www.eaglelabs.com
by e-mail: sales@eaglelabs.com
or by US Mail: Eagle Labs
10201-A Trademark Street
Rancho Cucamonga, CA 91730-5850 USA

Shipping within U.S. – FOB Rancho Cucamonga, California. Shipping and insurance will be added to the order invoice.

Payment Terms – Net 30 days

Returned product must be accompanied by an RGA (Returned Goods Authorization) Number. Only merchandise credit will be issued unless a shipping error occurs from the manufacturer (Eagle Labs). Returns are subject to a 10% re-stocking fee. Product returns will be accepted within one year from ship date.

All Eagle Products are “Made in the USA” and unconditionally guaranteed against manufacturing defects.

Eagle products are single use and provided sterile unless otherwise noted.

Quality Assurance



We at **Eagle Labs** guarantee that our products are manufactured with the highest quality standards. Eagle Labs complies with CFR 21 of The U.S. Food and Drug Administration. Our manufacturing facility is ISO 13485 and MDD 93/42/EEC Registered. Eagle Products also carry the European Union CE Mark Registration.

Commitment

Our commitment at Eagle Labs remains the same as it has been for over 25 years. We continually strive to offer high-quality and cost-effective Single-Use Medical Devices to the Ophthalmic Community. We proudly distribute our products for surgery throughout the world and support those healthcare professionals who have dedicated their time and efforts to improve quality of sight and life for everyone.

Table of Contents

Anesthesia Administration

Atkinson Tip Peribulbar Needles	5
Atkinson Tip Retrobulbar Needles	5
Sharp Point Retrobulbar Needles	5
Tri-Port Sub-Tenon Anesthesia Cannulae	5

Anterior Chamber Maintenance

Anterior Chamber Maintainer	5
Anterior Chamber Maintainer Kit	5
Anterior Chamber Maintainer, Pediatric	5

Capsule Polishing

Aspirating J-Hook Capsule Polisher	6
Capsule Polisher (Kratz)	6
Micro Loop Curette Polisher	6
Olive Tip Capsule Polisher	6
Silicone "Squeegee" Capsule Polishers	6

Capsulotomy

Angled Irrigating Cystotomes	8
Capsulorhexis Suction Cystotome (Dahan)	8
Curved Irrigating Cystotomes	8
Formed Irrigating Cystotomes	7, 8
Straight Irrigating Cystotomes	7, 8

Corneal Abrasion

Stromal Puncture Needle	8
-------------------------	---

Cortex Removal

Aspirating J-Hook Cannulae	9
Aspirating J-Hook Capsule Polisher	9
Aspirating C-Hook Cannulae	9
Aspirating J-Hook "Squeegee" Cannulae	9
Bimanual I/A Cortex Removal System	9
Cortex Aspiration Cannulae	9
Simcoe Cortex Extraction Cannulae	9

Drapes & Fluid Collection

Fluid Drainage Wick	11
Half Body Drapes	11
Medium Plastic Drapes	10
Micro Surgical Accessory Kit	11
Ophthalmic Surgery Kit	11
Phaco Drainage Bags	11
Small Plastic Drapes	10

Eye Protection

Bilateral Eye Shields, Clear	12
Corneal Light Shields	12
Eye Shields	12
Universal Eye Shields, Clear	12

Fluid Administration

Fluid Tubing Extension Sets	13
High Flow Filter	12
Fluid Administration Sets	13

Table of Contents (Continued)

Glaucoma

Anterior Chamber Cannulae	14
Glaucoma Revision Picks	14
Visco Canalostomy Cannulae	14

Hydrodissection & Hydrodelineation

Cortical Cleaving Hydrodissection Cannulae	15
Nucleus Hydrodissection Cannulae	14, 15
Nucleus Hydrodelineation Cannulae	16

IOL/Nucleus Manipulation

Intraocular Lens Manipulators	16
-------------------------------	----

Irrigation, Air Injection & Viscoelastic

Anterior Chamber Irrigation Cannulae	17
Anterior Chamber Washout Cannulae	17
Bishop Harmon Cannulae	17
FLOvisc Viscoelastic Cannulae	18
Infusion Cannulae	18
Stromal Hydration Cannulae	18
Tri-Step Viscoelastic Cannulae	18

Irrigation & Aspiration

Double I and A Cannulae	19
Double I and A Simcoe Cannulae	19

Knives & Blades **G O L D™** Premium Ophthalmic Surgical Knives

Blade Holders	34
Clear Corneal Knives	32
Coaxial MicroPhaco Knives	31
Corneal Teflon Block	31
Corneal Trephine Blades	31
Crescent Knives	33
Ear Blades	34
Lamellar Blades	35
Microsurgical Blades	34
Miniature Edged Blades	35
MVR Knives	32
MVR Knives, Angled	32
Needle Blades	35
Phaco Slit Knives	33
Scleral Blades	35
Short Stop IOL Inserter Knives	33
Spoon Knives	32
Stab Knives	33
Trapezoid Clear Corneal / Implant Knives	32

Nucleus Removal

Hydro-Expression Cannulae	20
Irrigating Vectis	20
Viscoelastic Expression Cannulae	20

Oculoplastic

Lacrimal Cannulae	.20
Lacrimal Intubation Sets (DCR)	.20

Refractive (LASIK)

PVA LASIK Eye Spears	.24
Pre-expanded LASIK Spears	.24
PVA LASIK Drainage Ring (Chayet)	.24
Irrigating LASIK Cannulae	.21, 22, 23
LASIK Drainage Ring	.24
LASIK Shields	.24
LASIK Sponge	.24
Small Refractive Drapes	.25

Speculums

Wire Eye Speculums	.25
--------------------	-----

Surgical Absorption Products

Cellulose Eye Sponges	.26
Diamond Knife Cleaning Block	.26
Eye Pledget	.26
Instrument Wipes	.26
Pre-Expanded LASIK Spears	.26
PVA Eye Spears	.26

Syringes, Handles & Adapters

Aspirating Handle	.26
Aspiration Syringes	.26
Irrigating Handle	.26
Male/Male Luer Adapter	.26
Spring Assisted Syringes	.26

Vitreoretinal


Infusion Cannulae	.30
High Viscosity Infusion Cannulae	.30
Micro-Vitreoretinal Cannulae	.29
Perfluorocarbon I/A Cannulae	.30
Retinal Dissectors and Micro Picks	.27
Silicone Oil Cannulae	.28
Silicone Oil Cannulae, Manual Delivery System	.27
Silicone Oil Removal Kit	.28
Subretinal Aspirating Cannulae	.29
Subretinal Extrusion Cannulae	.27
Subretinal Fluid Cannulae	.28
Vitreoretinal Injection Cannulae	.29
Vitreoretinal Viscodissectors	.30

Ordering Information	.1
-----------------------------	----

Product Index	.36-39
----------------------	--------

World Wide Distribution	.40
--------------------------------	-----

ANESTHESIA ADMINISTRATION




Atkinson Tip Retrobulbars

- 101-23** 23GA x 1 1/2" (38mm)
- 101-25** 25GA x 1 1/2" (38mm)
- 101-27** 27GA x 1 1/2" (38mm)

Atkinson Tip Peribulbars

- 102-23-32** 23GA x 1 1/4" (32mm)
- 102-25-32** 25GA x 1 1/4" (32mm)
- 102-27-32** 27GA x 1 1/4" (32mm)
- 102-23-22** 23GA x 7/8" (22mm)
- 102-25-22** 25GA x 7/8" (22mm)
- 102-27-22** 27GA x 7/8" (22mm)



Tri-Port Sub-Tenon Anesthesia Cannula

Can be used to administer antibiotics during Straubismus surgery (Mendez/Goldbaum)

2 - .4mm side ports, open-end port


- 100-19** 19GA x 1" angled (25mm)
- 100-19C** 19GA x 1" curved (25mm)
- 100-21** 21GA x 1" angled (25mm)
- 100-21C** 21GA x 1" curved (25mm)



Sharp Point Retrobulbars

- 103-23** 23GA x 1 1/2" (38mm)
- 103-25** 25GA x 1 1/2" (38mm)

ANTERIOR CHAMBER MAINTENANCE



3 per box

Anterior Chamber Maintainer Kit

- 170-20T** 20GA

(infusion cannula with self retaining serrated tip)
20GA x 3/16" (4.76mm)
cannula attached to 8" (20.3cm) flexible tubing
attached to an I.V. line.




3 per box

Anterior Chamber Maintainer

- 170-20** 20GA

(infusion cannula with self retaining serrated tip)
20GA x 3/16" (4.76mm) cannula
attached to 8" (20.3cm) flexible tubing



3 per box

Pediatric A.C. Maintainer

Used to maintain the chamber during pediatric retinal procedure (Mendez/Freeman)

- 171-20R** 20GA

3.5mm, 90° angled, 20GA cannula
with 40° beveled tip connected to 11" (28cm)
flexible tube (.8mm I.D., 1.6mm O.D.)

CAPSULE POLISHING

Silicone "Squeegee" Tip Capsule Polisher (Engels)

40 degree curve with extended silicone sleeve for polishing

128-27 27GA x 1" (25mm)





Silicone "Squeegee" Tip Capsule Polisher (Terry)

7mm bend to tip with flexible sleeve

125-23 23GA x 1" (25mm)

125-27 27GA x 1" (25mm)



Silicone "Squeegee" Tip Cortex Extractor and Polisher (Buedding)

Irrigating cortex extractor with flexible sleeve for the 12 o'clock position straight 23mm long

107-27J 27GA



Capsule Polisher (Kratz)

8mm bend, .30mm top port, rounded tip

125-21 21GA x 1" (25mm)

125-24 24GA x 1" (25mm)

Roughened end permits scraping and aspiration of cortical debris



Capsule Polisher (Kratz)

Roughened tip, angled 3mm

125-22 22GA x 1 1/2" (38mm)



Capsule Polisher (Lu)

124-23 23GA x 1" (25mm)

Curved, roughened olive tip - open port



Micro Loop Curette Polisher

125-25 25GA x 1 1/4" (32mm)

Irrigating posterior capsule polisher
curette loop is angled, posterior surface of loop is roughened, with open port



Olive Tip Capsule Polisher

10mm bend with roughened tip - open port

126-23 23GA x 1" (25mm)

126-25 25GA x 1" (25mm)



CAPSULOTOMY



Irrigating Cystotomes

13mm bend to tip (formed)

120-01 25GA x 3/4" (19mm)

120-02 27GA x 3/4" (19mm)



Irrigating Cystotomes/Push-Pull

13mm bend to tip (formed)

120-03 25GA x 3/4" (19mm)



Irrigating Cystotomes

Straight

120-21S 21GA x 1" (25mm)


120-22S 22GA x 1" (25mm)

120-23S 23GA x 1" (25mm)

120-25S 25GA x 5/8" (16mm)

120-27S 27GA x 1/2" (13mm)

120-30S 30GA x 3/4" (19mm)



Irrigating Cystotomes


15mm bend to tip (formed)

120-21F 21GA x 1" (25mm)

120-22F 22GA x 1" (25mm)

120-23F 23GA x 1" (25mm)

120-25F 25GA x 1" (25mm)



Irrigating Cystotomes


12mm bend to tip (formed)

120-22F-12 22GA x 3/4" (19mm)

120-23F-12 23GA x 3/4" (19mm)

120-25F-12 25GA x 5/8" (16mm)

120-27F-12 27GA x 1/2" (13mm)



10mm bend to tip (formed)

120-30F-10 30GA x 1/2" (13mm)



Irrigating Cystotomes

8mm bend to tip (formed)

120-25F-8 25GA x 5/8" (16mm)

120-27F-8 27GA x 1/2" (13mm)



Irrigating Cystotomes/Push-Pull (Blumenthal)

Straight

121-25S 25GA x 5/8" (16mm)


121-27S 27GA x 1/2" (13mm)



8mm bend to tip (formed)

121-25F-8 25GA x 5/8" (16mm)

121-27F-8 27GA x 1/2" (13mm)



12mm bend to tip (formed)

121-25F-12 25GA x 5/8" (16mm)

121-27F-12 27GA X 1/2" (13mm)

CAPSULOTOMY (continued)


Irrigating Cystotomes/Reverse Bend

	Straight
122-23S	23GA x 3/4" (19mm)
122-25S	25GA x 5/8" (16mm)
122-27S	27GA x 1/2" (13mm)
	12mm bend to tip (formed)
122-23RF	23GA x 3/4" (19mm)
122-25RF	25GA x 5/8" (16mm)
122-27RF	27GA x 1/2" (13mm)



Irrigating Cystotome (Joyce)

Curved for intercapsular/endocapsular techniques
220-25 25GA x 5/8" (16mm)



Irrigating Cystotome (Joyce)


Curved for intercapsular/endocapsular techniques
220-27 27GA x 1/2" (13mm)



Capsulorhexis Suction Cystotome (Dahan)

Curved design to create suction
 and hold capsular flap



221-25* 25GA x 1" (25mm)



*Used with extension tubing #169-30L-10 and syringe part #300-03 or #300-05

Irrigating Cystotomes/Angled



	Angled 45°, 10mm bend to tip
123-25	25GA x 5/8" (16mm)
123-27	27GA x 1/2" (13mm)



CORNEAL ABRASION

Stromal Puncture Needle


To achieve anterior stromal puncture of the epithelium
165-25 25GA x 1" (25mm)
 angled 30°, 19mm bend to tip
 sharp cutting tip



CORTEX REMOVAL



Irrigation



Aspiration




Bimanual I/A Cortex Removal System

192-23 Irrigation 21GA x 1" (25mm)
Two .4mm ports with female connector



Aspiration 23GA x 1" (25mm)
One .4mm port .5mm from tip
with male connector

Male and female connectors are compatible with I/A and Phaco systems. Allows cortex removal through paracentesis or Phaco incision while maintaining the chamber with irrigation.



Aspiration



Cortex Aspiration Cannula



191-21 21GA x 1" (25mm) .4mm port with male connector
191-23 23GA x 1" (25mm) .4mm port with male connector
191-25 25GA x 1" (25mm) .4mm port with male connector



Hook dimensions 2.5mm long x 2mm wide

Aspirating J-Hook Cannula

.4mm inside port
104-25 25GA x 1" (25mm)



Hook dimensions 2.75mm long x 2.5mm wide

Aspirating J-Hook Cannula/ Capsule Polisher (Thornton)

Open end port,
roughened on the end of hook
for capsule polishing

107-24 24GA x 1" (25mm)



Hook dimensions 2.0mm long x 2.0mm wide

Aspirating J-Hook Cannula

Angled 20mm long, straight 23mm long

105-22L, R, OR S 22GA

105-25L, R, OR S 25GA

105-27L, R, OR S 27GA

Left, right or straight hook



Hook dimensions 1.5mm long X 1.5mm wide

Aspirating Cortex Extractor

For aspirating cortex between the
10 and 2 o'clock positions

106-23S 23GA

106-25S 25GA

106-27L, R, OR S 27GA

106-30L, R, OR S 30GA

Left, right, or straight hook



Hook dimensions 1.5mm long X 2mm wide

Aspirating C-Hook Extractor

106-27SC 27GA



Simcoe Cortex Extractor Aspiration Cannula

Curved with .3mm inside port

108-23 23GA

108-25 25GA


DRAPES & FLUID COLLECTION

Small Plastic Drape w/ Adhesive Aperture

810-20

Overall size 16" x 16" (41cm x 41cm)
 Adhesive aperture area 5" x 5" (13cm x 13cm)
 Aperture circle 2.5" (6cm)

10 per box




Small Plastic Drape w/ Adhesive Incise

810-21

Overall size 16" x 16" (41cm x 41cm)
 Adhesive incise area 5" x 5" (13cm x 13cm)

10 per box




Refractive Drape w/ Split Aperture

810-22

Overall size 16" x 10.7" (40cm x 27cm)
 Adhesive aperture width 3" (7.6cm)
 Aperture oval width 1.5" (3.8cm)

10 per box




Medium Plastic Drape w/ Adhesive Aperture

810-30

Overall size 50" X 32" (127cm X 81cm)
 Adhesive aperture Area 6" X 5" (15cm X 13cm)
 Aperture oval 2" X 4" (5cm X 10cm)

5 per box




Medium Plastic Drape w/ Adhesive Incise

810-60

Overall size 48" X 50" (122cm X 127cm)
 Adhesive incise area 5" X 5" (13cm X 13cm)

5 per box




Medium Plastic Drape w/ Adhesive Incise & Fluid Collection Pouch


810-61

Overall size 48" X 51" (122cm X 130cm)
 Adhesive incise area 4" X 5" (10cm X 13cm)

5 per box




DRAPES & FLUID COLLECTION (continued)



5 per box

Non-Woven Half Body Drape w/ Adhesive Incise & Fluid Collection Pouch 810-54


Overall size 61" X 77" (155cm X 196cm)
Adhesive incise area 4" X 3" (10cm X 8cm)
w/ velcro cord holder



5 per box

Non-Woven Half Body Drape w/ Adhesive Aperture & Fluid Collection Pouch 810-64

Overall size 61" X 77" (155cm X 196cm)
Adhesive aperture area 4" X 3" (10cm X 8cm)
Aperture oval 1" X 2" (3cm X 5cm)
w/ velcro cord holder




20 per box

Phaco Drainage Bag - 400ml

3.875" x 8.25"; 400ml capacity with two drainage wicks and adhesive backing to adhere to surgical drape

700-03




20 per box

Drainage Bag - 80ml

2" x 8.25"; 80ml capacity with two drainage wicks and adhesive backing to adhere to surgical drape

700-02



20 per box

Microsurgical Accessory Kit

Cost Saving Convenience Kit includes one each:
400ml Drainage Bag with Two Wicks, Corneal Light Shield,
Instrument Wipe

700-06




20 per box

Fluid Drainage Wick

1/8" x 6.125" The biocompatible polymer is extremely hydrophilic and fast wicking, absorbing up to 25x its weight in fluids

700-04



20 per box

Ophthalmic Surgery Kit

Cost Saving Convenience Kit includes one each:
Wick, Corneal Light Shield, Instrument Wipe

700-07

EYE PROTECTION

Corneal Light Shield

8mm diameter. Protects the retina from light. Contains no cellulose and leaves no lint or debris. Upon moistening, becomes soft and porous, conforming to shape of cornea. Keeps cornea moist. Easy to remove.

700-01

20 per box



Universal Eye Shields, Clear

400-15

Clear polycarbonate lenses, vented, non-sterile

50 per box




Bilateral Eye Shields, Clear

400-25

Clear polycarbonate lenses, vented, w/foam lining & adjustable strap, non-sterile

25 per box




Eagle Eye Shields

400-01

Eye shield, non-sterile, left or right eye


25 per box



400-01S

Eye shield, sterile, left or right eye

10 per box



FLUID ADMINISTRATION

High Flow Filter .20 Micron

Male/female luer connections, used between syringe and cannula to filter out particulates and contaminants

201-20

50 per box



FLUID ADMINISTRATION (CONTINUED)



20 per box



20 per box


Fluid Tubing Sets

Extension sets with male/female luer connections provide fluid administration during ophthalmic procedures


169-30L-6 Extension set - 6" (15.2cm)
(I.D. - .8mm, O.D. -1.6mm)

169-30L-10 Extension set - 10" (25.4cm)
(I.D. - .8mm, O.D. -1.6mm)


169-30L-30 Extension set - 30" (76.2cm)
(I.D. - 3mm, O.D. - 4mm)




3 per box



3 per box



3 per box



3 per box

Fluid Administration Sets

Features an adjustable roller clamp for control of fluid flow

Single
200-01 Tubing length 86"

Double
200-02 Tubing length 96"

Single w/ .20 micron filter
200-21 Tubing length 86"


Double w/ .20 micron filter
200-22 Tubing length 96"

GLAUCOMA


Visco Canalostomy Cannula

177-33 5mm beveled tip, slightly angled

Micro-gauge cannula for viscoelastic injection during glaucoma surgery



5 per box




Glaucoma Revision Pick

60 Degree (Lederer)

177-27 27GA x 1/2" (12.5mm)

Angled 60 degree, 7mm bend to tip
tip has a sharp side cutting edge
(superior approach)



Glaucoma Revision Pick

30 Degree (Lederer)

178-27 27GA x 1/2" (12.5mm)

Angled 30 degree, 10mm bend to tip
tip has a sharp side cutting edge
(temporal or nasal approach)



Anterior Chamber Cannula (Mendez)

110-27TM 27GA x 10mm, round tip

Angled 45°, 3mm bend to tip

To reform the anterior chamber during glaucoma procedures



HYDRODISSECTION

Nucleus Hydrodissection Cannula

1" (25mm) long with flattened tip

113-23 23GA 7mm bend to tip

113-25 25GA 7mm bend to tip

113-27 27GA 7mm bend to tip

113-30 30GA 7mm bend to tip


113-25-8 25GA 8mm bend to tip

113-23-10 23GA 10mm bend to tip

113-25-10 25GA 10mm bend to tip

113-27-10 27GA 10mm bend to tip

113-30-10 30GA 10mm bend to tip



Hydrodissection Cannula (Tuebingen)

313-27 27GA x 1" (25mm), 75 degree curved, round tip

313-30 30GA x 1" (25mm), 75 degree curved, round tip




Nucleus Hydrodissection Cannula Vertical Flow

1" (25mm) long with flattened tip

213-25 25GA 7mm bend to tip

213-27 27GA 7mm bend to tip




Nucleus Hydrodissection Cannula (Nichamin)

1" (25mm) long with flattened tip, formed

113-25F-12 25GA 12mm bend to tip

113-27F-12 27GA 12mm bend to tip

113-30F-12 30GA 12mm bend to tip




HYDRODISSECTION (continued)



Nucleus Hydrodissection Cannula (Bobbitt)

Double angled cannula allows insertion into anterior capsule for hydrodissection

116-27 27GA x 1" (25mm)



Cortical Cleaving Hydrodissection Cannula

Angled 45 degree,
1" (25mm) long with round tip


216-30 30GA 12mm bend to tip



Cortical Cleaving Hydrodissection Cannula (Grabow)

Angled 45 degree,
1" (25mm) long with round tip

217-27 27GA 11mm bend to tip



Cortical Cleaving Hydrodissection Cannula

Angled 45 degree,
1" (25mm) long with flattened tip


218-27 27GA 11mm bend to tip



J-Hook Hydrodissection Cannula

For hydrodissection between the
10 - 2 o'clock positions

105-27FL 27GA
Hook dimensions 2.25mm long X 1.75mm wide,
flattened tip



Nucleus Hydrodissection Cannula (Johnson)

115-25 25GA x 1" (25mm)


115-27 27GA x 1" (25mm)

115-30 30GA x 1" (25mm)

10mm bend with beveled tip for ease of insertion
thru small incision.



horizontal flow



vertical flow

Nucleus Hydrodissection Cannula

118-25 25GA x 1" (25mm) curved (Helsinki) horizontal

118-25S 25GA x 1" (25mm) curved (Sauter) vertical

118-27 27GA x 1" (25mm) curved (Helsinki) horizontal

118-27S 27GA x 1" (25mm) curved (Sauter) vertical

HYDRODELINEATION

Nucleus Hydrodelineation Cannula (Johnson)

114-25 25GA x 1" (25mm)

114-27 27GA x 1" (25mm)

114-30 30GA x 1" (25mm)

With beveled tip for ease of insertion thru small incision

Nucleus Hydrodelineation Cannula

117-23 23GA


117-25 25GA

Parallel bevel of cannula allows hydrodynamic penetration of the cataract, separating the lenticular lamellae of the inner nucleus, without plugging the cannula lumen

Nucleus Hydrodelineator/Rotator

Formed 12mm from bend to tip. Tip is single bevel, angled 90°

119-25 25GA x 1" (25mm) curved



IOL / NUCLEUS MANIPULATION

Used for positioning intraocular lenses, handles have ridges for added control

IOL Manipulator (Jaffe/Maltzman)

127-01

Fits 0.2mm diameter holes or larger, angled 10mm



3 per box



IOL Manipulator (Sinskey)

127-03

Fits 0.2mm diameter holes or larger, straight



3 per box



IOL Manipulator (Sinskey)

127-04

Fits 0.2mm diameter holes or larger, angled 10mm



3 per box



Button Tip Manipulator

127-05

0.5mm button tip, angled 10mm

May also be used for iris retraction



3 per box



IOL Manipulator (Fenzl)

127-06

Fits 0.2mm diameter holes or larger, angled 10mm



3 per box



IRRIGATION, AIR INJECTION & VISCOELASTIC



Anterior Chamber Irrigator


1" (25mm) long with round tip, angled 30°

- 110-20** 20GA 7mm bend to tip
- 110-21** 21GA 7mm bend to tip
- 110-21-4** 21GA 4mm bend to tip
- 110-22** 22GA 7mm bend to tip
- 110-23** 23GA 7mm bend to tip
- 110-23-4** 23GA 4mm bend to tip
- 110-25-11** 25GA 11mm bend to tip
- 110-25** 25GA 7mm bend to tip
- 110-25-4** 25GA 4mm bend to tip
- 110-26** 26GA 7mm bend to tip
- 110-27-11** 27GA 11mm bend to tip
- 110-27-9** 27GA 9mm bend to tip
- 110-27** 27GA 7mm bend to tip
- 110-27-4** 27GA 4mm bend to tip
- 110-30-11** 30GA 11mm bend to tip
- 110-30** 30GA 7mm bend to tip
- 110-30-4** 30GA 4mm bend to tip



110-26^{1/2} 26GA 8mm bend to tip
1/2" (12mm) long with round tip, angled 25°

110-30SP 30GA 12mm bend to tip
9/16" (14mm) long with beveled tip, angled 45°



Anterior Chamber Washout Cannula


1" (25mm) long with flattened tip

- 111-16** 16GA
- 111-19** 19GA



Bishop Harmon Cannula

- 111-16BH** 16GA x 1" (25mm), 6mm bend to tip
- 111-19BH** 19GA x 1" (25mm), 8mm bend to tip



Anterior Chamber Irrigator, Straight

1" (25mm) long with round tip

- 112-20** 20GA
- 112-23** 23GA
- 112-25** 25GA
- 112-27** 27GA
- 112-30** 30GA



Anterior Chamber Irrigator, Curved

112-25C 25GA x 1" (25mm) curved with round tip


IRRIGATION, AIR INJECTION & VISCOELASTIC (CONTINUED)

FLOvisc™ CANNULAE

FLOvisc Cannulae are specifically designed for easier expression of viscoelastics. the internal insert eliminates any air pocket to reduce bubbling

Tri-Step FLOvisc Viscoelastic Cannula (DeCamp)


109-20-25 25GA tip
20GA cannula, stepped down to a 25GA cannula,
10mm bend to tip, angled 30°, 20GA cannula
Insert extends into the hub to further reduce bubbling



109-20-27 27GA tip
20GA cannula, stepped down to a 27GA cannula,
10mm bend to tip, angled 30°, 20GA cannula
Insert extends into the hub to further reduce bubbling

FLOvisc Viscoelastic Cannula (DeCamp)

109-22 27GA tip
22GA cannula, stepped down to a 27GA cannula,
10mm bend to tip, angled 30°, 27GA cannula
Insert extends into the hub to further reduce bubbling



109-27 27GA tip
20GA cannula, stepped down to a 27GA cannula,
10mm bend to tip, angled 30°



109-27-25 25GA tip
20GA cannula, stepped down to a 25GA cannula,
10mm bend to tip, angled 30°

Stromal Hydration Cannula

165-27 27GA x 1" (25mm), angled, 12mm bend
closed end flattened tip, inside port
Designed to irrigate the outer stroma for improved sealing of clear corneal incisions



Infusion Cannula (Murphy)

170-20M 20GA x 1" (25mm)
angled 45°, 7mm bend to tip,
2 EA. .4mm side ports,
45° beveled tip (up)





Infusion Cannula (Wright)

170-20W 20GA x 1" (25mm)
angled 45°, 7mm bend to tip,
2 EA. .4mm side ports,
45° beveled tip (down)




IRRIGATION & ASPIRATION



6 per box

I and A Cannula (Simcoe)



150-23 23GA x 23GA
Aspiration through 13" (33cm) flexible tubing with .4mm side port and irrigation through cannula



6 per box



Reverse I and A Cannula (Simcoe)

150-23R 23GA x 23GA
Aspiration through cannula with .4mm side port, irrigation through 13" (33cm) flexible tubing




Double I and A Cannula

130-25 23GA x 25GA x 1 1/2" (38mm)
angled 35°, 8mm from bend to tip




Double I and A Cannula

130-19 21GA x 19GA x 1 1/2" (38mm)
Irrigation/aspiration cannula
130-21 21GA x 21GA x 1 1/2" (38mm)
Irrigation/aspiration cannula




Double I and A Cannula

130-23 23GA x 23GA x 1 1/2" (38mm)
with 4mm angled aspiration tip



Double I and A Cannula

130-23L 23GA x 23GA x 1 1/2" (38mm)
Allows continuous irrigation/aspiration with iris retraction



Double I and A Cannula

130-23J 23GA x 23GA x 1 1/2" (38mm)
With J - hook aspiration tip

NUCLEUS REMOVAL

Viscoelastic Expression Cannula (Corydon)

109-25 25GA x 1" (25mm)

Curved design allows placement under the nucleus for viscoelastic expression during ECCE. Angled tip can be used to hook and remove nucleus.



Hydro-Expression Cannula (Corydon)

119-27 27GA x 1" (25mm)

Curved design allows placement under the nucleus for expression with balanced salt solution during ECCE. Angled tip can be used to hook and remove nucleus.



Irrigating Vectis

140-25-1 25GA, 5mm wide with 1 .3mm open port

140-25-2 25GA, 5mm wide with 2 .3mm open port

140-25-3 25GA, 5mm wide with 3 .3mm open port

***140-27-1** 27GA, 3mm wide with 1 .3mm open port

* For small incision surgery



OCULOPLASTIC

Lacrimal Cannula

Open end

160-21 21GA x 1 1/2" (38mm) curved

160-25 25GA x 1 1/2" (38mm) curved

160-26 26GA x 1 1/2" (38mm) curved




Closed end with side port

161-23 23GA x 1 1/2" (38mm) curved, .3mm port

161-25 25GA x 1 1/2" (38mm) curved, .3mm port

163-21 21GA x 1 1/2" (38mm) straight, .3mm port

163-25 25GA x 1 1/2" (38mm) straight, .3mm port



Lacrimal Cannula (Straight)

162-23 19GA x 1" (25mm) with 1/2" (13mm)

23GA tip

*Used to open up drainage of lacrimal system



Lacrimal Intubation Sets


12" (30.5cm) flexible tubing

O.D. .025" (.60mm), I.D. .012" (.30mm)

180-23 Probe 23GA x 7" (18cm)

180-27 Probe 27GA x 4 1/2" (11.5cm)


3 per box



180-23T Olive tip (Crawford)
(Same as 180-23 with olive tip)

180-27T Olive tip (Crawford)
(Same as 180-27 with olive tip)

3 per box



Lacrimal Intubation Sets, Curved (Guibor)

12" (30.5cm) flexible tubing.

O.D. .025" (.60mm), I.D. .012" (.30mm)

180-23C Probe 23GA x 7" (18cm)

3 per box




REFRACTIVE / LASIK




Linn LASIK Cannula

619-25 25GA x 1" (25mm) angled, 9mm bend
open end flattened tip, 2 ports




LASIK Irrigating Cannula (Banaji)

604-25 25GA x 1" (25mm)
Formed 12mm bend to tip,
6 ports, closed end



LASIK Irrigating Cannula (Guell)

605-25 25GA x 1" (25mm)
Formed 9mm bend to tip,
6 ports, open end




LASIK Irrigating Cannula (Burrato)

606-25 25GA x 1" (25mm)
Formed 9mm bend to tip,
2 ports, open end



LASIK Irrigating Cannula (Pettigrove)

607-25 25GA x 1" (25mm)
Formed 12mm bend to tip,
6 ports, open end




LASIK Irrigating Cannula (Zuckerman)

608-25 25GA x 1" (25mm)
angled 9mm bend to tip,
6 ports, open end



Irrigating LASIK Spatula (Slade)

613-26 26GA x 1" (25mm)
Formed 9mm bend to tip,
spatula tip



LASIK Flap Irrigator (Vidaaurri)

600-06 Dual 25GA multi-port Irrigation tips

This instrument features two angled parallel cannulae, each with four ports to simultaneously irrigate both sides of the LASIK flap and cleanse the stromal bed for flap repositioning.


REFRACTIVE / LASIK (CONTINUED)

Silicone Tip LASIK Flap Irrigator/Squeegee (Magruder)

602-23 23GA x 1" (25mm) angled, 7mm bend to tip

602-27 27GA x 1" (25mm) angled, 7mm bend to tip

Angled silicone tip allows for safe lifting of the flap, silicone sleeve acts as a "squeegee" to remove excess fluid and smooth down the flap to realign gutter




LASIK Flap Spatula/Irrigator (Hersh)

610-23 23GA x 1" (25mm) angled, 9mm bend to tip

610-25 25GA x 1" (25mm) angled, 9mm bend to tip


610-27 27GA x 1" (25mm) angled, 9mm bend to tip



611-23 23GA x 1" (25mm) angled, 4mm bend to tip

611-25 25GA x 1" (25mm) angled, 4mm bend to tip

611-27 27GA x 1" (25mm) angled, 4mm bend to tip



Non invasive spatula tip allows for safe lifting of the flap, curved angle for smoothing the flap during alignment, open end for irrigation

LASIK Irrigating Cannula (Lindstrom)

603-25 25GA x 1" (25mm) angled 7mm bend to tip

Features three each .3mm ports used to irrigate the anterior and posterior sides of the flap, smooth angle curvature allows for safe flap alignment




LASIK Fountain Cannula


600-21 25GA x 1" (25mm), right

600-22 25GA x 1" (25mm), left

Both fountain cannulae are 25GA with a single port and a 10mm curved angle to conform to the corneal shape. The tip is a smooth closed end design slightly tapered to allow insertion under the flap.



(Right)




(Left)



Quad Ported LASIK Irrigating Cannula

600-09 23GA x 1 1/2" (38mm), 4ea. .3mm ports

This instrument features four ports to irrigate the surgical interface and flush any debris collected during surgery



REFRACTIVE / LASIK (CONTINUED)



Nichamin **LASIK** Irrigating Cannula

- 600-23** 23GA x 1" (25mm) 12mm bend to tip
600-25 25GA x 1" (25mm) 12mm bend to tip
600-27 27GA x 1" (25mm) 12mm bend to tip

This instrument features a flattened tip with open-end irrigation, formed curvature used for elevating and smoothing the corneal flap



Nichamin **LASIK** Irrigating Cannula

- 601-23** 23GA x 1" (25mm) 6mm bend to tip
601-25 25GA x 1" (25mm) 6mm bend to tip
601-27 27GA x 1" (25mm) 6mm bend to tip

This instrument features a flattened tip with open-end irrigation, formed curvature used for elevating and smoothing the corneal flap



LASIK Aspiration Spoon


- 600-07** 25GA multi-port irrigation tip

This unique spoon-shaped instrument removes fluid and debris from beneath the flap



LASIK Irrigating Cannula

- 612-25** 25GA X 1" (25mm)
 Angled 45°, 12MM bend to tip



LASIK Irrigating Cannula (Soloman)

- 614-25** 25GA X 1" (25mm)
 Angled / formed spatula tip with 3 ports



LASIK Irrigating Cannula (Donnenfeld)

- 615-23** 23GA X 1" (25mm)
 Curved with a horizontal flattened tip



LASIK Irrigating Cannula

- 617-27** 27GA X 1" (25mm)
 Formed, 10mm bend to tip, flattened tip, 4 ports

REFRACTIVE / LASIK (CONTINUED)

LASIK Drains & Shields from Eagle Laboratories are made of the highest quality **LINT-FREE** PVA material to consistently ensure a particulate-free environment during LASIK Surgery. Drains are placed onto the orbital area after the microkeratome is used, before the flap is lifted. Used to control stromal bed hydration.

LASIK Drainage Ring (Chayet)

600-20 (12mm opening)

20 per box



LASIK Drainage Ring

600-30 Superior Hinge, (9.5mm opening)


20 per box



LASIK Sponge

600-40 (11mm opening)

20 per box



LASIK shields are placed over the posterior side of the flap and hinge, to moisten and protect against unwanted laser ablation.

LASIK Shields

600-08 Nasal Hinge, 7mm dia., (5mm edge)

20 per box



620-08 Nasal Hinge, 8mm dia., (6mm edge)

20 per box



620-09 Nasal Hinge, 9mm Dia., (7mm edge)

20 per box



620-10 Superior Hinge, Dual Edge, (4mm & 8mm edges)

20 per box



PVA Eye Spears

600-10

- **LINT-FREE** PVA material.
- Compressed spears are ideal for controlling and absorbing fluid in and around the orbital area during **LASIK**.

10 packs of 5




Pre-Expanded LASIK Spears


600-15

- Pre-expanded PVA eye spears offer a **LINT-FREE** absorption surface.
- Faster wicking than compressed spears.
- Pre-expanded spears are ideal for aligning and smoothing the corneal flap.

20 packs of 2




REFRACTIVE / LASIK (CONTINUED)



10 per box

Small Plastic Drape w/ Adhesive Aperture 810-20


Overall size	16" x 16" (41cm x 41cm)
Adhesive area	5" x 5" (13cm x 13cm)
Aperture, circle	2.5" (6cm)



10 per box

Small Plastic Drape w/ Adhesive Incise 810-21

Overall size	16" x 16" (41cm x 41cm)
Adhesive area	5" x 5" (13cm x 13cm)




10 per box

Refractive Drape w/ Split Aperture 810-22

Overall size	16" x 10.7" (40cm x 27cm)
Adhesive Aperture Width	3" (7.6cm)
Aperture oval width	1.5" (3.8cm)

SPECULUMS



3 per box

Wire Eye Speculum (Barraquer) 940-10 10mm



3 per box

Wire Eye Speculum (Barraquer) 940-15 15mm



3 per box

Wire Eye Speculum (Kratz / Barraquer) 940-20 15mm



3 per box

Wire Eye Speculum (Kratz / Barraquer) 940-29 15mm (temporal)

SURGICAL ABSORPTION PRODUCTS

Instrument Wipe

700-05

3" x 3"; medium density. Useful for cleaning delicate surgical instruments without leaving lint or debris.

20 per box



Diamond Knife Cleaning Block

700-15

Small pore construction absorbs fluids and removes debris from diamond blade when gentle pressure is applied.

10 per box



Cellulose Eye Sponges

500-01

High volume fast wicking action for absorption of blood and fluids within the operating field

10 packs of 10



PVA Eye Spears

600-10

- **LINT FREE** PVA material
- Ideal for corneal and refractive surgery

10 packs of 5




Pre-Expanded LASIK Spears

600-15

- Pre-expanded PVA eye spears offer a **LINT-FREE** absorption surface.
- Faster wicking than compressed spears.
- Pre-expanded spears are ideal for aligning and smoothing the corneal flap.

20 packs of 2



Eye Pledget

500-05

2mm x 6mm

LINT FREE PVA Sponge

10 packs of 10



SYRINGES, HANDLES & ADAPTERS

Spring Assisted Aspirating Syringe (Pallin)

300-03

3ml aspirating syringe with Luer Lok

300-05

5ml aspirating syringe with Luer Lok

10 per box



Aspirating Syringe

300-03SP

3ml aspirating syringe with Luer Lok

300-05SP

5ml aspirating syringe with Luer Lok

10 per box



Male/Male Luer Adapter

191-01

Male luer fittings on both ends, used for connecting Phaco and I/A tubing to instruments.

50 per box



Irrigating Handle

191-02

3" (7.5cm)

Luer fittings both ends

10 per box



Aspirating Handle

191-03

3" (7.5cm)

Luer fittings both ends, with aspiration port

10 per box



VITREORETINAL



Membrane Dissector (Kingham)

169-21 21GA x 1" (25mm)

Sharp tip with cutting edge angled 45°



Vitreoretinal Micro Pick (Glaser)

169-25 20GA x 1" (25mm) with 25GA x 3/16" (4.76mm) tip

Sharp tip with cutting edge angled 50°



Membrane Micro Pick (Eaton)

169-30 20GA x 1" (25mm) with 30GA x 3/16" (4.76mm) tip

Sharp tip with cutting edge angled 50°



Epiretinal Membrane Dissector (Wagner)

168-21 21GA x 1" (25mm)

Sharp tip with cutting edge angled 45°



Vitreoretinal Micro Pick


168-25 20GA x 1" (25mm) with 25GA x 3/16" (4.76mm) tip

Sharp tip with cutting edge angled 50°



Hard Tapered Extrusion Cannula (Tim Mendez)

176-25 20GA x 1" with 25GA x 1mm tip



Silicone Oil Cannula, Manual Delivery System

With 6" (15cm) high pressure extension tube pre-attached to a 10 cc ring control syringe

(Mendez/Freeman)

167-20TS 20GA x 8mm

167-19TS 19GA x 8mm

(Wagner)


167-16TS 16GA x 8mm

VITREORETINAL (CONTINUED)

Silicone Oil Cannula (Wagner)

167-16T 16GA x 8mm with 6" (15cm)
High pressure extension tube

3 per box




Silicone Oil Cannula (Peyman)

167-16 16GA x 8mm

167-19 19GA x 8mm

For injection or aspiration of high viscosity fluids in the repair of retinal detachments



Silicone Oil Removal Kit (T. Mendez) (Sork)*

80" PVC Tubing with connectors and 19GA straight, Cannula fits both the Alcon and Storz cassettes

164-19 19GA

3 per box

* Allows the safe removal of silicone oil in 10 minutes or less




Subretinal Fluid Cannula 6mm (Grizzard)

20GA cannula w/ silicone sleeve extended 6mm beyond end, overall length with tubing 37mm.

166-25 straight

166-25C curved

Soft silicone tip ensures safe entry through a macular tear or hole for aspiration of subretinal fluids.




Subretinal Fluid Cannula 3mm (Grizzard)

20GA cannula w/ silicone sleeve extended 3mm beyond end, overall length with tubing 34mm.


166-20 straight

166-21 curved

Soft silicone tip ensures safe entry through a macular tear or hole for aspiration of subretinal fluids.



VITREORETINAL (CONTINUED)



Vitreoretinal Injection Cannula Straight (Vander)

172-30 20GA x 1" (25mm)
With 30GA x 5mm round tip



Vitreoretinal Injection Cannula Angled (Vander)

172-30C 20GA x 1" (25mm)
With 30GA round tip



Micro Vitreoretinal Cannula Angled (Glaser)

171-32 20GA x 1" (25mm)
With 32GA tip, 4mm long,
angled 65° 1mm from tip



Subretinal Injection Cannula Straight (Glaser)

173-32 20GA x 1" (25mm)
With 32GA x 3mm, 35° beveled tip




Subretinal Injection Cannula Curved (Glaser)

173-32C 20GA x 1" (25mm)
With 32GA x 3mm, 35° beveled tip



Subretinal Aspiration Cannula (Landers)

174-27 20GA x 1" (25mm)
With 27GA, 6.5mm tip, angled
2.5mm AT 60° from 22.5° beveled tip.



Vitreoretinal Aspiration Cannula (Landers)

174-28 20GA x 1" (25mm)
With 27GA x 5/32" (4mm) round tip



Subretinal Aspiration Cannula (Lambert)

173-34 20GA x 1" (25mm)
With 32GA x 5/16" (8mm) curved
beveled tip



Vitreoretinal Injection Cannula (Mendez)


175-32 20GA x 1 1/4" (32mm)
With 32GA (7mm) curved round tip

VITREORETINAL (CONTINUED)

Perfluorocarbon Co-axial I/A Cannula

Used for perfluorocarbon liquid exchange during retinal procedures.


- 149-19** 25GA X 1 1/2" (38mm)
With open end irrigation port
19GA outer dia. with .4mm aspiration port



Perfluorocarbon I/A Cannula

Used for perfluorocarbon liquid exchange during retinal procedures.

- 149-23** 23GA X 1 1/2" (38mm)
With open end irrigation port
25GA aspiration cannula




- 149-25** 25GA X 1 1/2" (38mm)
With open end irrigation port
25GA aspiration cannula

Infusion Cannula

- 175-25** 20GA X 2.5mm, beveled tip
175-40 20GA X 4mm, beveled tip
175-60 20GA X 6mm, beveled tip

High Viscosity Infusion Cannula

- 175-14** 20GA X 4mm, beveled tip
175-16 20GA X 6mm, beveled tip



3 per box





Vitreoretinal Viscodissector (Glaser)

- 169-30-90** 20GA x 1" (25mm) with 30GA x 3/16" (4.76mm) tip 90°
169-30-70 20GA x 1" (25mm) with 30GA x 3/16" (4.76mm) tip 70°
169-30-45 20GA x 1" (25mm) with 30GA x 3/16" (4.76mm) tip 45°

30GA cannula angled 90°, 70° OR 45° with 1mm flattened tip. Cannula is attached to a Luer lock high pressure extension tube (.03" I.D. x 8" long).

3 per box



KNIVES AND BLADES



1 per box

Corneal Trephine Blades

707-00	7.00mm
707-25	7.25mm
707-50	7.50mm
707-75	7.75mm
708-00	8.00mm
708-25	8.25mm
708-50	8.50mm
708-75	8.75mm



1 per box

Corneal Teflon Block

701-01




6 per box



GOLD™
Premium Ophthalmic Surgical Knives
BRAND

Coaxial MicroPhaco Knives

EG-1122	2.2 mm, angled, Full handle, double bevel
----------------	--



6 per box




GOLD™
Premium Ophthalmic Surgical Knives
BRAND

Sub 2.0mm Phaco Slit Knives

For bi-manual Phaco
and micro-incision cataract procedures

EG-1014	1.4mm, angled, full handle, bevel up
EG-1015	1.5mm, angled, full handle, bevel up
EG-1016	1.6mm, angled, full handle, bevel up
EG-1017	1.7mm, angled, full handle, bevel up
EG-1018	1.8mm, angled, full handle, bevel up
EG-1020	2.0mm, angled, full handle, bevel up



6 per box



GOLD™
Premium Ophthalmic Surgical Knives
BRAND

Double Bevel Sideport Knives

To create a self-sealing incision
during bi-manual Phaco procedures

EG-5010	1.0mm, angled, full handle, double bevel
EG-5012	1.2mm, angled, full handle, double bevel
EG-5014	1.4mm, angled, full handle, double bevel

KNIVES AND BLADES (CONTINUED)

Spoon Knives

For scleral tunnel incision

EG-6821 3.0 mm, angled,
full handle, double bevel



6 per box



GOLD™
Premium Ophthalmic Surgical Knives
BRAND

MVR Knives

For vitreoretinal infusion port
or paracentesis incision

EG-5555 19GA, full handle

EG-5560 20GA, full handle

EG-5563 23GA, full handle



6 per box




GOLD™
Premium Ophthalmic Surgical Knives
BRAND

MVR Knives Angled

For vitreoretinal infusion port
or paracentesis incision

EG-5660 20GA, full handle, angled

EG-5663 23GA, full handle, angled



6 per box



GOLD™
Premium Ophthalmic Surgical Knives
BRAND

Clear Corneal Knives

EG-1127 2.75 mm, angled,
full handle, double bevel

EG-1130 3.0 mm, angled,
full handle, double bevel

EG-1132 3.2 mm, angled,
full handle, double bevel



6 per box



GOLD™
Premium Ophthalmic Surgical Knives
BRAND

Clear Corneal Knives, Parallel Sides

EG-1327 2.75 mm angled,
full handle, double bevel

EG-1330 3.0 mm angled,
full handle, double bevel

EG-1332 3.2 mm angled,
full handle, double bevel



6 per box




GOLD™
Premium Ophthalmic Surgical Knives
BRAND

Trapezoid Clear Corneal / Implant Knives

EG-2832 2.85 mm - 3.2 mm angled,
full handle, bevel up

EG-3032 3.0 mm - 3.2 mm angled,
full handle, bevel up

Dual-purpose design provides the ability to
accomplish clear corneal and implant incisions
with one knife.



6 per box



GOLD™
Premium Ophthalmic Surgical Knives
BRAND

KNIVES AND BLADES

6 per box



GOLD™
Premium Ophthalmic Surgical Knives BRAND

Stab Knives

For side port or scleral incision

- EG-1500** 15°, full handle
- EG-2250** 22.5°, full handle
- EG-3000** 30°, full handle
- EG-4500** 45°, full handle

6 per box




GOLD™
Premium Ophthalmic Surgical Knives BRAND

Phaco Slit Knives

For clear corneal or scleral incision

- EG-5510** 2.5 mm, angled, full handle, bevel up
- EG-5568** 2.65 mm, angled, full handle, bevel up
- EG-4738** 2.75 mm, angled, full handle, bevel up
- EG-4740** 2.85 mm, angled, full handle, bevel up
- EG-4741** 2.85 mm, straight, full handle
- EG-5570** 3.0 mm, angled, full handle, bevel up
- EG-5571** 3.0 mm, straight, full handle
- EG-5572** 3.2 mm, angled, full handle, bevel up
- EG-5573** 3.2 mm, straight, full handle
- EG-5580** 3.5 mm, angled, full handle, bevel up

6 per box



GOLD™
Premium Ophthalmic Surgical Knives BRAND

Short Stop IOL Inserter Knives

For IOL enlargement

- EG-9336** 3.5 mm, angled, full handle, bevel up
- EG-9340** 4.0 mm, angled, full handle, bevel up
- EG-9350** 5.2 mm, angled, full handle, bevel up

6 per box



GOLD™
Premium Ophthalmic Surgical Knives BRAND

Crescent Knives

For Scleral tunnel incision


- EG-9402** 2.0 mm, angled, full handle, bevel up

KNIVES AND BLADES (CONTINUED)

Microsurgical Blades (Flat Stock)

EL-7511	Orange	1.5mm, 15°	blade and tip
EL-7513	Blue	3.0mm, 15°	blade and tip
EL-7515	Green	5.0mm, 15°	blade and tip
EL-7522	Grey	4.0mm, 22.5°	blade and tip
EL-7530	Black	3.5mm, 30°	blade and tip
EL-7545	Tan	3.0mm, 45°	blade and tip

6 units per box




Precise microsurgical blades for controlled depth incisions

Microsurgical Blades Full Handle (Flat Stock)

EL-7512	Orange	1.5mm, 15°
EL-7514	Blue	3.0mm, 15°
EL-7516	Green	5.0mm, 15°
EL-7523	Grey	4.0mm, 22.5°
EL-7531	Black	3.5mm, 30°
EL-7546	Tan	3.0mm, 45°

6 units per box




The safety and convenience offered only by one-piece, single use microsurgical knives, unsurpassed depth control and cutting angles

Ear Blades

EL-7100	Myringotomy
EL-7120	Myringotomy spear cut
EL-7121	Myringotomy lance cut
EL-7200	Spoon 60°

6 units per box



For ENT procedures

Blade Holders

EL-4312	7.5cm universal handle
EL-4310	10cm universal handle

1 unit per box



KNIVES AND BLADES (CONTINUED)



Versatile blades useful in pediatric surgery, capsulotomy, and paracenteses

Needle Blades

EL-5910	3.0mm
EL-5920	5.0mm
EL-5930	7.0mm
EL-5940	1.5mm
6 units per box	



Flat, polished undersurfaces make ideal undermining blades

Lamellar Blades

EL-6600	60°
6 units per box	



Full Handle Scleral Blades

EL-6823H	
6 units per box	



Scleral Blades

EL-5700	
6 units per box	



A variety of shapes and edges for precise microsurgical procedures

Miniature Edged Blades

EL-6200	Chisel tip
EL-6300	Pointed tip
EL-6400	Round tip
EL-6500	Lance tip
EL-6700	Curved tip
EL-6900	Sharp all around
20 units per box - 6200, 6400, 6500, & 6700	
6 units per box - 6300 & 6900	

Product Index

Product Number	Description	Box Quantity	Page	Product Number	Description	Box Quantity	Page
100-19	Tri-Port Sub-Tenon Cannula	10	5	110-25	Anterior Chamber Irrigator	10	17
100-19C	Tri-Port Sub-Tenon Cannula	10	5	110-26	Anterior Chamber Irrigator	10	17
100-21	Tri-Port Sub-Tenon Cannula	10	5	110-26-1/2	Anterior Chamber Irrigator	10	17
100-21C	Tri-Port Sub-Tenon Cannula	10	5	110-27	Anterior Chamber Irrigator	10	17
101-23	Retrobulbar Atkinson Tip	10	5	110-27-9	Anterior Chamber Irrigator	10	17
101-25	Retrobulbar Atkinson Tip	10	5	110-27TM	Anterior Chamber Cannula (Mendez)	10	14
101-27	Retrobulbar Atkinson Tip	10	5	110-30	Anterior Chamber Irrigator	10	17
102-23-22	Peribulbar Atkinson Tip	10	5	110-25-4	Anterior Chamber Irrigator	10	17
102-25-22	Peribulbar Atkinson Tip	10	5	110-27-4	Anterior Chamber Irrigator	10	17
102-27-22	Peribulbar Atkinson Tip	10	5	110-30-4	Anterior Chamber Irrigator	10	17
102-23-32	Peribulbar Atkinson Tip	10	5	110-25-11	Anterior Chamber Irrigator	10	17
102-25-32	Peribulbar Atkinson Tip	10	5	110-27-11	Anterior Chamber Irrigator	10	17
102-27-32	Peribulbar Atkinson Tip	10	5	110-30-11	Anterior Chamber Irrigator	10	17
103-23	Retrobulbar Sharp Point	10	5	110-30SP	Anterior Chamber Irrigator	10	17
103-25	Retrobulbar Sharp Point	10	5	111-16	Washout Cannula	10	17
104-25	Aspirating J-Hook Cannula	10	9	111-19	Washout Cannula	10	17
105-22L	Aspirating J-Hook Cannula	10	9	111-16BH	Bishop Harmon Cannula	10	17
105-22R	Aspirating J-Hook Cannula	10	9	111-19BH	Bishop Harmon Cannula	10	17
105-22S	Aspirating J-Hook Cannula	10	9	112-20	Anterior Chamber Irrigator	10	17
105-25L	Aspirating J-Hook Cannula	10	9	112-23	Anterior Chamber Irrigator	10	17
105-25R	Aspirating J-Hook Cannula	10	9	112-25	Anterior Chamber Irrigator	10	17
105-25S	Aspirating J-Hook Cannula	10	9	112-25C	Anterior Chamber Irrigator, Curved	10	17
105-27L	Aspirating J-Hook Cannula	10	9	112-27	Anterior Chamber Irrigator	10	17
105-27R	Aspirating J-Hook Cannula	10	9	112-30	Anterior Chamber Irrigator	10	17
105-27S	Aspirating J-Hook Cannula	10	9	113-23	Nucleus Hydrodissection Cannula	10	14
105-27FL	J-Hook Hydrodissection Cannula	10	15	113-25	Nucleus Hydrodissection Cannula	10	14
106-23S	Aspirating Cortex Extractor	10	9	113-27	Nucleus Hydrodissection Cannula	10	14
106-25S	Aspirating Cortex Extractor	10	9	113-30	Nucleus Hydrodissection Cannula	10	14
106-27L	Aspirating Cortex Extractor	10	9	113-25-8	Nucleus Hydrodissection Cannula	10	14
106-27R	Aspirating Cortex Extractor	10	9	113-23-10	Nucleus Hydrodissection Cannula	10	14
106-27S	Aspirating Cortex Extractor	10	9	113-25-10	Nucleus Hydrodissection Cannula	10	14
106-27SC	Aspirating Cortex Extractor	10	9	113-27-10	Nucleus Hydrodissection Cannula	10	14
106-30L	Aspirating Cortex Extractor	10	9	113-30-10	Nucleus Hydrodissection Cannula	10	14
106-30R	Aspirating Cortex Extractor	10	9	113-25F-12	Nucleus Hydrodissection Cannula	10	14
106-30S	Aspirating Cortex Extractor	10	9	113-27F-12	Nucleus Hydrodissection Cannula	10	14
107-24	Aspirating J-Hook Cannula	10	9	113-30F-12	Nucleus Hydrodissection Cannula	10	14
107-27J	Squeegee Cortex Extractor	10	6	114-25	Nucleus Hydrodelination Cannula	10	16
108-23	Simcoe Cortex Extractor	10	9	114-27	Nucleus Hydrodelination Cannula	10	16
108-25	Simcoe Cortex Extractor	10	9	114-30	Nucleus Hydrodelination Cannula	10	16
109-20-25	Tri-Step Viscoelastic Cannula	10	18	115-25	Nucleus Hydrodissection Cannula	10	15
109-20-27	Tri-Step Viscoelastic Cannula	10	18	115-27	Nucleus Hydrodissection Cannula	10	15
109-22	FLOvisc Viscoelastic Cannula	10	18	115-30	Nucleus Hydrodissection Cannula	10	15
109-25	Viscoelastic Expression Cannula	10	20	116-27	Nucleus Hydrodissection Cannula	10	15
109-27	FLOvisc Viscoelastic Cannula	10	18	117-23	Nucleus Hydrodelination Cannula	10	16
109-27-25	FLOvisc Viscoelastic Cannula	10	18	117-25	Nucleus Hydrodelination Cannula	10	16
110-20	Anterior Chamber Irrigator	10	17	118-25	Nucleus Hydrodissection Cannula	10	15
110-21	Anterior Chamber Irrigator	10	17	118-25S	Nucleus Hydrodissection Cannula	10	15
110-21-4	Anterior Chamber Irrigator	10	17	118-27	Nucleus Hydrodissection Cannula	10	15
110-22	Anterior Chamber Irrigator	10	17	118-27S	Nucleus Hydrodissection Cannula	10	15
110-23	Anterior Chamber Irrigator	10	17	119-25	Nucleus Hydrodelineator/Rotator	10	16
110-23-4	Anterior Chamber Irrigator	10	17	119-27	Hydro-Expression Cannula	10	20

Product Index (Continued)

Product Number	Description	Box Quantity	Page	Product Number	Description	Box Quantity	Page
120-01	Irrigating Cystotome	10	7	130-23	Double I/A Cannula	10	19
120-02	Irrigating Cystotome	10	7	130-23J	Double I/A Cannula	10	19
120-03	Irrigating Cystotome	10	7	130-23L	Double I/A Cannula	10	19
120-21S	Irrigating Cystotome	10	7	130-25	Double I/A Cannula	10	19
120-22S	Irrigating Cystotome	10	7	140-25-1	Irrigating Vectis Loop	10	20
120-23S	Irrigating Cystotome	10	7	140-25-2	Irrigating Vectis Loop	10	20
120-25S	Irrigating Cystotome	10	7	140-25-3	Irrigating Vectis Loop	10	20
120-27S	Irrigating Cystotome	10	7	140-27-1	Irrigating Vectis Loop	10	20
120-30S	Irrigating Cystotome	10	7	149-19	Perfluorocarbon Co-axial I/A Cannula	10	30
120-21F	Irrigating Cystotome	10	7	149-23	Perfluorocarbon I/A Cannula	10	30
120-22F	Irrigating Cystotome	10	7	149-25	Perfluorocarbon I/A Cannula	10	30
120-23F	Irrigating Cystotome	10	7	150-23	Simcoe I/A Cannula	6	19
120-25F	Irrigating Cystotome	10	7	150-23R	Simcoe I/A Cannula Reverse	6	19
120-25F-8	Irrigating Cystotome	10	7	160-21	Lacrimonal Cannula	10	20
120-27F-8	Irrigating Cystotome	10	7	160-25	Lacrimonal Cannula	10	20
120-22F-12	Irrigating Cystotome	10	7	160-26	Lacrimonal Cannula	10	20
120-23F-12	Irrigating Cystotome	10	7	161-23	Lacrimonal Cannula	10	20
120-25F-12	Irrigating Cystotome	10	7	161-25	Lacrimonal Cannula	10	20
120-27F-12	Irrigating Cystotome	10	7	162-23	Lacrimonal Cannula	10	20
120-30F-10	Irrigating Cystotome	10	7	163-21	Lacrimonal Cannula	10	20
121-25S	Irrigating Cystotome/Push Pull	10	7	163-25	Lacrimonal Cannula	10	20
121-27S	Irrigating Cystotome/Push Pull	10	7	164-19	Silicone Oil Removal Kit	3	28
121-25F-8	Irrigating Cystotome/Push Pull	10	7	165-25	Stromal Puncture Needle	10	8
121-27F-8	Irrigating Cystotome/Push Pull	10	7	165-27	Stromal Hydration Cannula	10	18
121-25F-12	Irrigating Cystotome/Push Pull	10	7	166-20	Subretinal Fluid Cannula, 3mm	10	28
121-27F-12	Irrigating Cystotome/Push Pull	10	7	166-21	Subretinal Fluid Cannula, 3mm	10	28
122-23S	Irrigating Cystotome/Reverse Bend	10	8	166-25	Subretinal Fluid Cannula, 6mm	10	28
122-25S	Irrigating Cystotome/Reverse Bend	10	8	166-25C	Subretinal Fluid Cannula, 6mm	10	28
122-27S	Irrigating Cystotome/Reverse Bend	10	8	167-16	Silicone Oil Cannula	10	28
122-23RF	Irrigating Cystotome/Reverse Bend	10	8	167-16T	Silicone Oil Cannula	3	28
122-25RF	Irrigating Cystotome/Reverse Bend	10	8	167-16TS	Silicone Oil Cannula	3	27
122-27RF	Irrigating Cystotome/Reverse Bend	10	8	167-19	Silicone Oil Cannula	10	27
123-25	Irrigating Cystotome/Angled	10	8	167-19TS	Silicone Oil Cannula	3	27
123-27	Irrigating Cystotome/Angled	10	8	167-20TS	Silicone Oil Cannula	3	27
124-23	Capsule Polisher (Lu)	10	6	168-21	Membrane Dissector	10	27
125-21	Capsule Polisher	10	6	168-25	Vitreoretinal Micro Pick	10	27
125-22	Capsule Polisher	10	6	169-21	Membrane Dissector	10	27
125-23	Squeegee Capsule Polisher	10	6	169-25	Vitreoretinal Micro Pick	10	27
125-24	Capsule Polisher	10	6	169-30	Membrane Micro Pick	10	27
125-25	Micro Loop Polisher	10	6	169-30-90	Vitreoretinal Viscodissector	3	30
125-27	Squeegee Capsule Polisher	10	6	169-30-70	Vitreoretinal Viscodissector	3	30
126-23	Olive Tip Polisher	10	6	169-30-45	Vitreoretinal Viscodissector	3	30
126-25	Olive Tip Polisher	10	6	169-30L-6	Fluid Tubing Extension Set	20	13
127-01	IOL Manipulator	3	16	169-30L-10	Fluid Tubing Extension Set	20	13
127-03	IOL Manipulator	3	16	169-30L-30	Fluid Tubing Extension Set	20	13
127-04	IOL Manipulator	3	16	170-20	Anterior Chamber Maintainer	3	5
127-05	Button Tip Manipulator	3	16	170-20M	Infusion Cannula	10	18
127-06	IOL Manipulator	3	16	170-20T	Anterior Chamber Maintainer Kit	3	5
128-27	Squeegee Capsule Polisher	10	6	170-20W	Infusion Cannula	10	18
130-19	Double I/A Cannula	10	19	171-20R	Anterior Chamber Maintainer	3	5
130-21	Double I/A Cannula	10	19	171-32	Micro Vitreoretinal Cannula	10	29

Product Index (Continued)

Product Number	Description	Box Quantity	Page	Product Number	Description	Box Quantity	Page
172-30	Vitreoretinal Injection Cannula	10	29	400-25	Bilateral Eye Shields, Clear	25	12
172-30C	Vitreoretinal Injection Cannula	10	29	500-01	Eye Sponges	10 packs of 10	26
173-32	Subretinal Injection Cannula	10	29	500-05	Eye Pledget	10 packs of 10	26
173-32C	Subretinal Injection Cannula	10	29	600-06	LASIK Flap Irrigator	10	21
173-34	Subretinal Injection Cannula (Lambert)	10	29	600-07	LASIK Aspiration Cannula	10	23
174-27	Subretinal Injection Cannula	10	29	600-08	LASIK Shield	20	24
174-28	Vitreoretinal Aspiration Cannula (Landers)	10	29	600-09	LASIK Irrigation Cannula	10	22
175-14	High Viscosity Infusion Cannula	3	30	600-10	PVA Eye Spears	10 packs of 5	24,26
175-16	High Viscosity Infusion Cannula	3	30	600-15	PVA Pre-Expanded LASIK Spears	20 packs of 2	24,26
175-25	Infusion Cannula 2.5mm	3	30	600-20	LASIK Drainage Ring, (Chayet)	20	24
175-32	Vitreoretinal Injection Cannula	10	29	600-21	LASIK Fountain Cannula,Right	10	22
175-40	Infusion Cannula 4mm	3	30	600-22	LASIK Fountain Cannula,Left	10	22
175-60	Infusion Cannula 6mm	3	30	600-23	LASIK Irrigation Cannula	10	23
176-25	Extrusion Cannula	10	27	600-25	LASIK Irrigation Cannula	10	23
177-27	Glaucoma Revision Pick	10	14	600-27	LASIK Irrigation Cannula	10	23
178-27	Glaucoma Revision Pick	10	14	600-30	LASIK Drainage Ring (9.5mm opening)	20	24
177-33	Visco Canalostomy Cannula	5	14	600-40	LASIK Sponge	20	24
180-23	Lacrimal Intubation Set	3	20	601-23	LASIK Irrigation Cannula	10	23
180-23C	Lacrimal Intubation Set	3	20	601-25	LASIK Irrigation Cannula	10	23
180-23T	Lacrimal Intubation Set	3	20	601-27	LASIK Irrigation Cannula	10	23
180-27	Lacrimal Intubation Set	3	20	602-23	Silicone Tip LASIK Irrigator	10	22
180-27T	Lacrimal Intubation Set	3	20	602-27	Silicone Tip LASIK Irrigator	10	22
191-01	Male/Male Luer Adapter	50	26	603-25	LASIK Irrigation Cannula	10	22
191-02	Irrigating Handle	10	26	604-25	LASIK Irrigating Cannula	10	21
191-03	Aspirating Handle	10	26	605-25	LASIK Irrigating Cannula	10	21
191-21	Cortex Aspiration Cannula	10	9	606-25	LASIK Irrigating Cannula	10	21
191-23	Cortex Aspiration Cannula	10	9	607-25	LASIK Irrigating Cannula	10	21
191-25	Cortex Aspiration Cannula	10	9	608-25	LASIK Irrigating Cannula	10	21
192-23	Bimanual I/A Cortex Removal System	10	9	610-23	LASIK Flap Spatula/Irrigator	10	22
200-01	Fluid Administration Sets Single	3	13	610-25	LASIK Flap Spatula/Irrigator	10	22
200-02	Fluid Administration Sets Double	3	13	610-27	LASIK Flap Spatula/Irrigator	10	22
200-21	Fluid Administration Sets Single Port w/Filter	3	13	611-23	LASIK Flap Spatula/Irrigator	10	22
200-22	Fluid Administration Sets Double Port w/Filter	3	13	611-25	LASIK Flap Spatula/Irrigator	10	22
201-20	High Flow Filter .20 Micron	50	12	611-27	LASIK Flap Spatula/Irrigator	10	22
213-25	Hydrodissection Cannula	10	14	612-25	LASIK Irrigating Cannula	10	23
213-27	Hydrodissection Cannula	10	14	613-26	Irrigating LASIK Spatula	10	21
216-30	Hydrodissection Cannula	10	15	614-25	LASIK Irrigating Cannula	10	23
217-27	Hydrodissection Cannula	10	15	615-23	LASIK Irrigating Cannula	10	23
218-27	Hydrodissection Cannula	10	15	617-27	LASIK Irrigating Cannula	10	23
220-25	Irrigating Cystotome	10	8	619-25	Linn LASIK Cannula	10	21
220-27	Irrigating Cystotome	10	8	620-08	LASIK Shield, 8mm Dia., (6mm edge)	20	24
221-25	Suction Cystotome	10	8	620-09	LASIK Shield, 9mm Dia., (7mm edge)	20	24
300-03	Spring Assisted Aspirating Syringe, 3ml	10	26	620-10	LASIK Shield, Dual Edge 4mm & 8mm	20	24
300-05	Spring Assisted Aspirating Syringe, 5ml	10	26	700-01	Corneal Light Shield	20	12
300-03SP	Aspirating Syringe, 3ml	10	26	700-02	Phaco Drainage Bag 80ml	20	11
300-05SP	Aspirating Syringe, 5ml	10	26	700-03	Phaco Drainage Bag 400ml	20	11
313-27	Hydrodissection Cannula (Tuebingen)	10	14	700-04	Fluid Drainage Wick	20	11
313-30	Hydrodissection Cannula (Tuebingen)	10	14	700-05	Instrument Wipe	20	26
400-01	Eye Shield, Non-Sterile	25	12	700-06	Microsurgical Accessory Kit	20	11
400-01S	Eye Shield	10	12	700-07	Ophthalmic Surgery Kit	20	11
400-15	Universal Eye Shield, Clear (non-sterile)	50	12	700-15	Diamond Knife Cleaning Block	10	26

Product Index (Continued)

Product Number	Description	Box Quantity	Page	Product Number	Description	Box Quantity	Page
701-01	Corneal Teflon Block	1	31	EG-5570	Phaco Slit Knife	6	33
707-00	Corneal Trephine Blade	1	31	EG-5571	Phaco Slit Knife	6	33
707-25	Corneal Trephine Blade	1	31	EG-5572	Phaco Slit Knife	6	33
707-50	Corneal Trephine Blade	1	31	EG-5573	Phaco Slit Knife	6	33
707-75	Corneal Trephine Blade	1	31	EG-5580	Phaco Slit Knife	6	33
708-00	Corneal Trephine Blade	1	31	EG-5660	MVR Knife, Angled	6	32
708-25	Corneal Trephine Blade	1	31	EG-5663	MVR Knife, Angled	6	32
708-50	Corneal Trephine Blade	1	31	EG-6821	Spoon Knife	6	32
708-75	Corneal Trephine Blade	1	31	EG-9336	Short Stop IOL Inserter Knife	6	33
810-20	Small Plastic Drape w/Adhesive Aperature	10	10,25	EG-9340	Short Stop IOL Inserter Knife	6	33
810-21	Small Plastic Drape w/Adhesive Incise	10	10,25	EG-9350	Short Stop IOL Inserter Knife	6	33
810-22	Refractive Drape w/ Split Aperture	10	10,25	EG-9402	Crescent Knife	6	33
810-30	Medium Plastic Drape w/Aperture	5	10	EL-4310	10cm Universal Handle, Non-Sterile	1	34
810-54	Half Body Drape w/Incise	5	11	EL-4312	7.5cm Universal Handle, Non-Sterile	1	34
810-60	Medium Plastic Drape w/Incise	5	10	EL-5700	Scleral Blade, CM	6	35
810-61	Medium Plastic Drape w/Incise	5	10	EL-5910	3.0mm Needle Blade, CM	6	35
810-64	Half Body Drape w/Aperture	5	11	EL-5920	5.0mm Needle Blade, CM	6	35
940-10	Wire Eye Speculum	3	25	EL-5930	7.0mm Needle Blade, CM	6	35
940-15	Wire Eye Speculum	3	25	EL-5940	1.5mm Needle Blade, CM	6	35
940-20	Wire Eye Speculum	3	25	EL-6200	Chisel Tip Blade, CM	20	35
940-29	Wire Eye Speculum	3	25	EL-6300	Pointed Tip Blade, CM	6	35
EG-1014	Sub 2.0mm Phaco Slit Knife	6	31	EL-6400	Round Tip Blade, CM	20	35
EG-1015	Sub 2.0mm Phaco Slit Knife	6	31	EL-6500	Lance Tip Blade, CM	20	35
EG-1016	Sub 2.0mm Phaco Slit Knife	6	31	EL-6600	60 Degree Lamellar Blade, CM	6	35
EG-1017	Sub 2.0mm Phaco Slit Knife	6	31	EL-6700	Curved Tip Blade, CM	20	35
EG-1018	Sub 2.0mm Phaco Slit Knife	6	31	EL-6823H	Scleral Blade, FH	6	35
EG-1020	Sub 2.0mm Phaco Slit Knife	6	31	EL-6900	Sharp Tip All Around, CM	6	35
EG-1122	Coaxial MicroPhaco Knife	6	31	EL-7100	Myringotomy Blade, CM	6	34
EG-1127	Clear Corneal Knife	6	32	EL-7120	Myringotomy Blade, CM	6	34
EG-1130	Clear Corneal Knife	6	32	EL-7121	Myringotomy Blade, CM	6	34
EG-1132	Clear Corneal Knife	6	32	EL-7200	Myringotomy Blade, CM	6	34
EG-1327	Clear Corneal Knife,Parallel Sides	6	32	EL-7511	1.5mm x 15 Degree Blade, CM	6	34
EG-1330	Clear Corneal Knife,Parallel Sides	6	32	EL-7512	1.5mm x 15 Degree Blade, FH	6	34
EG-1332	Clear Corneal Knife,Parallel Sides	6	32	EL-7513	3.0mm x 15 Degree Blade, CM	6	34
EG-1500	Stab Knife	6	33	EL-7514	3.0mm x 15 Degree Blade, FH	6	34
EG-2250	Stab Knife	6	33	EL-7515	5.0mm x 15 Degree Blade, CM	6	34
EG-2832	Trapezoid Clear Corneal/ Implant Knife	6	32	EL-7516	5.0mm x 15 Degree Blade, FH	6	34
EG-3000	Stab Knife	6	33	EL-7522	4.0mm x 22.5 Degree Blade, CM	6	34
EG-3032	Trapezoid Clear Corneal/ Implant Knife	6	32	EL-7523	4.0mm x 22.5 Degree Blade, FH	6	34
EG-4500	Stab Knife	6	33	EL-7530	3.5mm x 30 Degree Blade, CM	6	34
EG-4738	Phaco Slit Knife	6	33	EL-7531	3.5mm x 30 Degree Blade, FH	6	34
EG-4740	Phaco Slit Knife	6	33	EL-7545	3.0mm x 45 Degree Blade, CM	6	34
EG-4741	Phaco Slit Knife	6	33	EL-7546	3.0mm x 45 Degree Blade, FH	6	34
EG-5010	Double Bevel Sideport Knife	6	31				
EG-5012	Double Bevel Sideport Knife	6	31				
EG-5014	Double Bevel Sideport Knife	6	31				
EG-5510	Phaco Slit Knife	6	33				
EG-5555	MVR Knife	6	32				
EG-5560	MVR Knife	6	32				
EG-5563	MVR Knife	6	32				
EG-5568	Phaco Slit Knife	6	33				

- **Argentina**
Rosinov Surgical
Instruments (GSJ S.A.)
Fax 11-48553668 (ext. 222)
Tel. 11-48553668
www.silviorosinov.com
- **Australia**
Southern Cross
Hospital Supplies
Fax 02-9869-8716
Tel. 02-9890-5088
- **Austria / Eastern Europe**
Askin & Co.
Fax 1-9798846
Tel. 1-9798844
www.askin.co.at
- **Belgium**
DeCo-Vision
Fax 0032-9-245-47-73
Tel. 0032-9-245-47-72
- **Brazil**
Mediphacos
Fax 31-2102-2212
Tel. 31-2102-2211
www.mediphacos.com
- **Chile**
Oftomed
Fax 56-2-378-3445
Tel. 56-2-376-9210
www.oftomed.cl
- **Costa Rica**
Insuned Inc.
Fax 2-296-7424
Tel. 2-296-7422
- **Denmark**
Lindcon Optical Group A/S
Fax 86-99-37-77
Tel. 86-99-22-77
www.lindcon.dk
- **Egypt**
Medi-Tec
Fax 2-5-728-717
Tel. 2-5-717-844
- **France**
Carl Zeiss Meditec SAS
Fax 5-46-44-85-60
Tel. 5-46-44-85-50
www.meditec.zeiss.com/iol
- **Germany**
M+C
Fax 6182-929292
Tel. 6182-92920
www.m-und-c.de
- **Greece**
Pharmex S.A.
Fax 1-5144279
Tel. 1-5150781
www.pharmex.gr
- **Indonesia**
Universal Star
Multilink (USM)
Fax 62-21-29001419
Tel. 62-21-29001517-19
usm.jhalim@gmail.com
- **Italy**
Alfa Intes Casoria
Fax 081-584-2526
Tel. 081-759-6622
www.alfaintes.it
- **Japan**
Nihon Surge
Fax 6-6190-7867
Tel. 6-6190-7865
- **Latvia**
Latvijas-Amerikas
Acu Centrs
Fax 7-272128
Tel. 7-272257
www.laac.lv
- **Mexico**
Oftalmica Internacional
www.oftalmica.com

Tijuana office
Fax 664-634-04-41
Tel. 01800-02-66065
marcela@oftalmica.com

Mexico City office
Fax 55-5574-94-96
Tel. 01800-01-52800
oisamex@hotmail.com
- **Norway**
Medistim Norge AS
Fax 47-23-03-52-51
Tel. 47-23-03-52-50
www.medistim.com
- **Perú**
Droguería ASG Medical
Fax 1-446-5259
Tel. 1-446-5259
- **Poland**
Consultronix S.A.
Fax 12-290-22-23
Tel. 12-290-22-25
www.consultronix.pl
- **Saudi Arabia**
Arab Supply Est.
Fax 1-463-1516
Tel. 1-279-1526
www.arabsupply.com.sa

FGC Medical
Fax 1-2191888
Tel. 1-2190888
www.fgcmedical.com
- **South Africa**
National Medical Supplies
Fax 012-3865576
Tel. 012-3865698
www.nationalmedical.co.za
- **South Korea**
Sung Woo Medical Co., Ltd.
Fax 2-785-6513
Tel. 2-783-4052
- **Spain**
Bloss Group
Fax 93-308-49-50
Tel. 93-303-65-65
www.blossgroup.com
- **Sweden**
Argus Medicinteknik AB
Fax 0480-49-13-80
Tel. 0480-49-12-80
www.argusmedicinteknik.se
- **Switzerland**
medilas ag
Fax 44-747-40-05
Tel. 44-747-40-00
www.medilas.ch
- **Turkey**
Care Optik Market Ltd.
Fax 212-233-40-99
Tel. 212-296-76-46
www.optikmarket.com
- **U.K. / Ireland**
Dixey Instruments
Fax 1-604-882488
Tel. 1-604-882480
www.dixeyinstruments.com
- **United Kingdom**
M.I.S.S.
Fax 153-646-0762
Tel. 153-674-1139
- **Uruguay**
Alejandro Besón
Fax 598-2336-0567
Tel. 598-2336-0567



10201-A Trademark Street
Rancho Cucamonga
California 91730-5850 USA

Toll Free: 800 782-6534

Ph: 909 481-0011

Fax: 909 481-4481

e-mail: goeagle@eaglelabs.com

www.eaglelabs.com

© Eagle Labs 2013

